
Product Information

Product SKU:	CAB10516	Gene ID:	320538	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	190kDa	Conjugate:	Unconjugated
Calculated MW:	130kDa/131kDa/141kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in DNA replication-independent chromatin assembly. Predicted to be located in nucleoplasm. Predicted to be active in nucleus. Is expressed in central nervous system and retina. Orthologous to human UBN2 (ubinnuclein 2).

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: mKIAA2030; 2900060J04Rik; 6030408G03Rik; D130059P03Rik; UBN2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse UBN2 (NP_796159.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.