

NEPRO Rabbit Polyclonal Antibody



CAB10520

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

63kDa

Calculated MW:

63kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

212547

Uniprot

Synonyms:

NEPRO; C3orf17; NET17; bc027231

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 263-564 of mouse NEPRO (NP_666084.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - NEPRO Polyclonal Antibody (CAB10520)

