

CAB1059

Product Information

Product SKU: CAB1059	Gene ID: 64223	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 37kDa	Conjugate: Unconjugated
Calculated MW: 36kDa	Isotype: IgG

Immunogen Information

Background: Enables protein serine/threonine kinase activator activity. Involved in TORC1 signaling; positive regulation of TOR signaling; and regulation of actin cytoskeleton organization. Part of TORC1 complex and TORC2 complex.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: GBL; LST8; POP3; WAT1; GbetaL; [KD Validated] MLST8

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human MLST8 (NP_001186103.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.