

VKORC1L1 Rat Polyclonal Antibody



CAB10600

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

24KDa

Calculated MW:

19kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rat

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

69568

Uniprot

Synonyms:

VKORC1L1


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of mouse VKORC1L1 (NP_081397.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - VKORC1L1 Rat pAb (CAB10600)