

RELL2 Rabbit Polyclonal Antibody



CAB10620

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

32kDa

Applications:

IHC

Reactivity:

Human

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

285613

Uniprot

Q8NC24

Synonyms:

RELL2; C5orf16

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 44-303 of human RELL2 (NP_776189.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
