

MT1B Rabbit Polyclonal Antibody



CAB10633

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

6kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

The protein encoded by this gene binds heavy metals and protects against toxicity from heavy metal ions. This gene is found in a cluster of similar genes on chromosome 16.

Immunogen information

Gene ID:

4490

Uniprot

P07438

Synonyms:

MT1B; MT-1B; MT-IB; MT1; MT1Q; MTP

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-61 of human MT1B (NP_005938.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
