

# SPC25 Rabbit Polyclonal Antibody



CAB10653

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

26kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

This gene encodes a protein that may be involved in kinetochore-microtubule interaction and spindle checkpoint activity.

## Immunogen information

### Gene ID:

57405

### Uniprot

Q9HBM1

### Synonyms:

SPC25; AD024; SPBC25; hSpc25

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human SPC25 (NP\_065726.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---