

AHI1 Rabbit Polyclonal Antibody



CAB10678

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

70kDa/120kDa/137kDa

Applications:

IHC

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified.

Immunogen information

Gene ID:

54806

Uniprot

Q8N157

Synonyms:

AHI1; AHI-1; JBTS3; ORF1; dJ71N10.1; jouberin

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human AHI1 (NP_001128304.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
