

# RNASE12 Rabbit Polyclonal Antibody



CAB1072

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

17kDa

### Applications:

IHC

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

493901

#### Uniprot

Q5GAN4

#### Synonyms:

RNASE12; HEL-S-85p; RAI1

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human RNASE12 (NP\_001019993.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---