

CAB10773

Product Information

Product SKU:	CAB10773	Gene ID:	25861	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	97kDa	Isotype:	IgG

Immunogen Information

Background:	This gene is thought to function in the organization and stabilization of stereocilia elongation and actin cytoskeletal assembly, based on studies of the related mouse gene. Mutations in this gene have been associated with autosomal recessive non-syndromic deafness and Usher Syndrome. Alternative splicing results in multiple transcript variants.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	WI; CIP98; USH2D; DFNB31; PDZD7B; WHRN
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 708-907 of human WHRN (NP_056219.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.