

CAB1079

Product Information

Product SKU:	CAB1079	Gene ID:	81788	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

Background: Enables ATP binding activity; magnesium ion binding activity; and protein serine/threonine kinase activity. Involved in several processes, including cellular response to glucose starvation; negative regulation of apoptotic process; and protein phosphorylation. Predicted to be active in cytoplasm. Implicated in anencephaly.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ANPH2; SNARK; NUAK2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 479-628 of human NUAK2 (NP_112214.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.