

CAB11352

Product Information

Product SKU:	CAB11352	Gene ID:	1956	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	175kDa	Conjugate:	Unconjugated
Calculated MW:	134kDa	Isotype:	IgG

Immunogen Information

Background:	The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor, thus inducing receptor dimerization and tyrosine autophosphorylation leading to cell proliferation. Mutations in this gene are associated with lung cancer. EGFR is a component of the cytokine storm which contributes to a severe form of Coronavirus Disease 2019 (COVID-19) resulting from infection with severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2).
Recommended Dilution:	WB,1:500 - 1:1000 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cellsIF/ICC,1:50 - 1:200
Synonyms:	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; EGFR
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1021-1210 of human EGFR (NP_005219.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.