

CAB11481

Product Information

Product SKU:	CAB11481	Gene ID:	2167	Size:	20uL, 100uL
Clone No:	ARC0616	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	15kDa	Conjugate:	Unconjugated
Calculated MW:	15kDa	Isotype:	IgG

Immunogen Information

Background: FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	aP2; ALBP; AFABP; A-FABP; HEL-S-104; FABP4
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FABP4 (P15090).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.