

FABP4 Rabbit Monoclonal Antibody



CAB11481

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

15KDa

Calculated MW:

15kDa

Applications:

WB

Reactivity:

Mouse, Rat

Protein Background

FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism. [provided by RefSeq, Jul 2008]

Immunogen information

Gene ID:

2167

Uniprot

P15090

Synonyms:

A-FABP; AFABP; ALBP; HEL-S-104; aP2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A synthesized peptide derived from human FABP4

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - FABP4 Rabbit mAb (CAB11481)