

CAB11488

Product Information

Product SKU:	CAB11488	Gene ID:	2247	Size:	20uL, 100uL
Clone No:	ARC0618	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	18kDa/22kDa/31kDa	Conjugate:	Unconjugated
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

Background:	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members bind heparin and possess broad mitogenic and angiogenic activities. This protein has been implicated in diverse biological processes, such as limb and nervous system development, wound healing, and tumor growth. The mRNA for this gene contains multiple polyadenylation sites, and is alternatively translated from non-AUG (CUG) and AUG initiation codons, resulting in five different isoforms with distinct properties. The CUG-initiated isoforms are localized in the nucleus and are responsible for the intracrine effect, whereas, the AUG-initiated form is mostly cytosolic and is responsible for the paracrine and autocrine effects of this FGF.
Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:100 - 1:500 IF/ICC,1:50 - 1:200
Synonyms:	BFGF; FGFB; FGF-2; HBGF-2; F2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human FGF2 (NP_001997.5).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.