

CAB11680

Product Information

Product SKU:	CAB11680	Gene ID:	563	Size:	20uL, 100uL
Clone No:	ARC0673	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Involved in cell adhesion and detection of chemical stimulus involved in sensory perception of bitter taste. Located in extracellular space.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	ZAG; ZA2G; Zinc-α2-glycoprotein (ZAG/AZGP1)
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Zinc-α2-glycoprotein (ZAG/AZGP1) (P25311).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.