

ALKBH5 Rabbit Polyclonal Antibody



CAB11684

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

40-52KDa

Calculated MW:

43kDa/44kDa/52kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54890

Uniprot

Q6P6C2

Synonyms:

ALKBH5; ABH5; OFOXD; OFOXD1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 245-394 of human ALKBH5 (NP_060228.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

