

ENHO Rabbit Polyclonal Antibody



CAB11733

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

7kDa

Applications:

IHC

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

375704

Uniprot

Q6UWT2

Synonyms:

ENHO; C9orf165; UNQ470; adropin


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-76 of human ENHO (NP_940975.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry - ENHO Polyclonal Antibody (CAB11733)