

**CAB11735**

---

**Product Information**

<b>Product SKU:</b>	CAB11735	<b>Gene ID:</b>	1399	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	37kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	34kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes a protein kinase containing SH2 and SH3 (src homology) domains which has been shown to activate the RAS and JUN kinase signaling pathways and transform fibroblasts in a RAS-dependent fashion. It is a substrate of the BCR-ABL tyrosine kinase, plays a role in fibroblast transformation by BCR-ABL, and may be oncogenic.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** CRKL

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human CRKL (NP\_005198.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.