

Product Information

Product SKU:	CAB11759	Gene ID:	214105	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	84kDa	Isotype:	IgG

Immunogen Information

Background:	Enables DNA-binding transcription activator activity, RNA polymerase II-specific and DNA-binding transcription repressor activity, RNA polymerase II-specific. Involved in negative regulation of transcription by RNA polymerase II; positive regulation of transcription by RNA polymerase II; and proacrosomal vesicle fusion. Located in chromocenter and nucleus. Is expressed in several structures, including brain; genitourinary system; gut; heart; and lung. Orthologous to human SOX30 (SRY-box transcription factor 30).
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SOX30
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of mouse SOX30 (NP_775560.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.