

ADGRB2 Rabbit Polyclonal Antibody



CAB11957

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

169kDa/171kDa/172kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a seven-span transmembrane protein that is thought to be a member of the secretin receptor family. The encoded protein is a brain-specific inhibitor of angiogenesis. The mature peptide may be further cleaved into additional products (PMID:20367554). Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

576

Uniprot

O60241

Synonyms:

ADGRB2; BAI2; AGRB2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 700-800 of human ADGRB2 (NP_001281264.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
