

# Gsdmc2 Rabbit Polyclonal Antibody



CAB12002

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

54kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

331063

#### Uniprot

Q2KHK6

#### Synonyms:

GSDMC; MLZE; AI987692; 9130023H10; gasdermin-C; Gsdmc2

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 367-467 of mouse Gsdmc2 (NP\_001161746.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---