

CAB1204

Product Information

Product SKU:	CAB1204	Gene ID:	158062	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	18kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable small molecule binding activity. Predicted to be involved in single fertilization. Predicted to be located in extracellular region.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	LCN5; hLcn5; UNQ643; LCN6
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 25-120 of human LCN6 (NP_945184.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.