

Zswim8 Rabbit Polyclonal Antibody



CAB12056

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

197kDa

Calculated MW:

70kDa/88kDa/193kDa/196kDa/197kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

268721

Uniprot

Synonyms:

ZSWIM8; KIAA0913; Zswim8


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1600-1700 of mouse Zswim8 (NP_082272.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - Zswim8 Polyclonal Antibody (CAB12056)