

# SPON2 Rabbit Polyclonal Antibody



CAB12077

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

36kDa

### Calculated MW:

35kDa

### Applications:

WB IHC

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

10417

#### Uniprot

Q9BUD6

#### Synonyms:

SPON2; DIL-1; DIL1; M-SPONDIN; MINDIN; spondin-2

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 27-150 of human SPON2 (NP\_036577.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - SPON2 Rabbit pAb (CAB12077)

