

Product Information

Product SKU:	CAB12077	Gene ID:	10417	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable antigen binding activity; lipopolysaccharide binding activity; and metal ion binding activity. Predicted to be involved in cell adhesion. Predicted to act upstream of or within several processes, including defense response to other organism; opsonization; and positive regulation of cytokine production. Located in extracellular exosome.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: DIL1; DIL-1; MINDIN; M-SPONDIN; SPON2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 27-150 of human SPON2 (NP_036577.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.