

RNMT Rabbit Polyclonal Antibody



CAB12103

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

55kDa

Calculated MW:

41kDa/46kDa/53kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC

1:100 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

67897

Uniprot

Synonyms:

RNMT; CMT1; CMT1c; MET; Met; RG7MT1; cm1p; hCMT1; hMet; 2610002P10Rik; AI848273; Rg7mt1; mKIAA0398


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of mouse RNMT (NP_080716.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - RNMT Polyclonal Antibody (CAB12103)