

SYCP1 Rabbit Polyclonal Antibody

CAB12139



Product Information

Product SKU:	CAB12139	Gene ID:	6847	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	120kDa	Conjugate:	Unconjugated
Calculated MW:	114kDa	Isotype:	IgG

Immunogen Information

Background:	Enables double-stranded DNA binding activity. Involved in protein homotetramerization. Predicted to be located in synaptonemal complex. Predicted to be active in central element; male germ cell nucleus; and transverse filament.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CT8; SCP1; SCP-1; HOM-TES-14; SYCP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human SYCP1 (NP_003167.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.