

CAB12192

Product Information

Product SKU:	CAB12192	Gene ID:	79814	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable agmatinase activity. Predicted to be involved in putrescine biosynthetic process from arginine, using agmatinase. Predicted to be located in mitochondrion.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	AGMAT
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-352 of human AGMAT (NP_079034.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.