

CAB12229

Product Information

Product SKU:	CAB12229	Gene ID:	10787	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	125kDa	Conjugate:	Unconjugated
Calculated MW:	129kDa	Isotype:	IgG

Immunogen Information

Background: Contributes to small GTPase binding activity. Involved in Rac protein signal transduction; positive regulation of Arp2/3 complex-mediated actin nucleation; and positive regulation of lamellipodium assembly. Located in extracellular exosome and focal adhesion. Part of SCAR complex.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: HEM2; NAP1; NAP125; p125Nap1; NCKAP1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human NCKAP1 (NP_038464.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.