

KIAA1524 Rabbit Polyclonal Antibody

CAB12267



Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB12267 | Gene ID: | 57650 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|--------|------------|--------------|
| Observed MW: | 102kDa | Conjugate: | Unconjugated |
| Calculated MW: | 102kDa | Isotype: | IgG |

Immunogen Information

Background: Enables protein homodimerization activity. Predicted to act upstream of or within positive regulation of neural precursor cell proliferation and spermatogenesis. Located in cytosol and plasma membrane.

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:1000 |
| Synonyms: | p90; NOCIVA; KIAA1524 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 681-905 of human KIAA1524 (NP_065941.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |