

CAB12271

Product Information

Product SKU:	CAB12271	Gene ID:	2079	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	12kDa	Conjugate:	Unconjugated
Calculated MW:	12kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables methyl-CpG binding activity. Predicted to be involved in pyrimidine nucleoside metabolic process. Located in midbody and nucleus. Part of methylosome.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	DROER; ERH
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-104 of human ERH (NP_004441.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.