

Product Information

Product SKU:	CAB1229	Gene ID:	5562	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	64kDa	Conjugate:	Unconjugated
Calculated MW:	64kDa	Isotype:	IgG

Immunogen Information

Background:	The protein encoded by this gene belongs to the ser/thr protein kinase family. It is the catalytic subunit of the 5'-prime-AMP-activated protein kinase (AMPK). AMPK is a cellular energy sensor conserved in all eukaryotic cells. The kinase activity of AMPK is activated by the stimuli that increase the cellular AMP/ATP ratio. AMPK regulates the activities of a number of key metabolic enzymes through phosphorylation. It protects cells from stresses that cause ATP depletion by switching off ATP-consuming biosynthetic pathways. Alternatively spliced transcript variants encoding distinct isoforms have been observed.
Recommended Dilution:	WB,1:1000 - 1:5000 IHC-P,1:100 - 1:500 IF/ICC,1:50 - 1:200
Synonyms:	AMPK; AMPKa1; AMPK alpha 1; AMPKα1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 381-520 of human AMPKα1 (NP_006242.5).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.