

PI 3 Kinase Class 3 Rabbit Monoclonal Antibody



CAB12295

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

100KDa

Calculated MW:

100kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5289

Uniprot

Q8NEB9

Synonyms:

VPS34; Vps34; hVps34


Immunogen:

A synthesized peptide derived from human PI 3 Kinase Class 3

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - PI 3 Kinase Class 3 Rabbit mAb (CAB12295)