

NDP Rabbit Polyclonal Antibody



CAB12470

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

15kDa

Applications:

IHC

Reactivity:

Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a secreted protein with a cystein-knot motif that activates the Wnt/beta-catenin pathway. The protein forms disulfide-linked oligomers in the extracellular matrix. Mutations in this gene result in Norrie disease and X-linked exudative vitreoretinopathy.

Immunogen information

Gene ID:

4693

Uniprot

Q00604

Synonyms:

NDP; EVR2; FEVR; ND; norrin


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human NDP (NP_000257.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded rat brain using NDP antibody (CAB12470) at dilution of 1:100 (40x lens).