TRPC1 Rabbit Polyclonal Antibody





Product Information

CAB12525

Product SKU: CAB12525 **Gene ID**: 7220 **Size**: 20uL, 100uL

Clone No: - Host Species: Rabbit Reactivity: Human, Mouse

Additional Information

Observed MW: 90kDa **Conjugate:** Unconjugated

Calculated MW: 91kDa Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a membrane protein that can form a non-selective channel

permeable to calcium and other cations. The encoded protein appears to be induced to form channels by a receptor tyrosine kinase-activated phosphatidylinositol second messenger system and also by

depletion of intracellular calcium stores. Two transcript variants encoding different isoforms have been

found for this gene.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: TRP1; HTRP-1; TRPC1

Purifcation Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPC1

(NP_001238774.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.