

# EVC Rabbit Polyclonal Antibody



CAB12634

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

111kDa

### Applications:

WB

### Reactivity:

Human

## Protein Background

This gene encodes a protein containing a leucine zipper and a transmembrane domain. This gene has been implicated in both Ellis-van Creveld syndrome (EvC) and Weyers acrorenal dysostosis.

## Immunogen information

### Gene ID:

2121

### Uniprot

P57679

### Synonyms:

EVC; DWF-1; EVC1; EVCL

## Antibody Information

### Recommended dilutions:

WB 1:1000 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 743-992 of human EVC (NP\_714928.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---