

ZC3H12A Rabbit Polyclonal Antibody



CAB12668

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

65kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:1000 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

80149

Uniprot

Q5D1E8

Synonyms:

ZC3H12A; MCPIP; MCPIP-1; MCPIP1; Reg1; dJ423B22.1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZC3H12A (NP_079355.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
