

C1QTNF4 Rabbit Polyclonal Antibody



CAB12730

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

35kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:1000 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

114900

Uniprot

Q9BXJ3

Synonyms:

C1QTNF4; CTRP4; ZACRP4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 250-329 of human C1QTNF4 (NP_114115.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
