

IFNA8 Rabbit Polyclonal Antibody



CAB12746

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

21kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

3445

Uniprot

P32881

Synonyms:

IFNA8; IFN-alphaB

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human IFNA8 (NP_002161.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
