

Product Information

Product SKU:	CAB12752	Gene ID:	51552	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	24kDa	Conjugate:	Unconjugated
Calculated MW:	24kDa	Isotype:	IgG

Immunogen Information

Background:	RAB14 belongs to the large RAB family of low molecular mass GTPases that are involved in intracellular membrane trafficking. These proteins act as molecular switches that flip between an inactive GDP-bound state and an active GTP-bound state in which they recruit downstream effector proteins onto membranes (Junutula et al., 2004 [PubMed 15004230]).
Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	FBP; RAB-14; RAB14
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 39-140 of human RAB14 (NP_057406.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.