

ABI Rabbit Polyclonal Antibody



CAB12786

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

73kDa

Applications:

WB

Reactivity:

Zea mays

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

542523

Uniprot

Synonyms:

ABI; VP1

Immunogen:

A synthetic peptide of Zea mays VP1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
