

# EIN3 Rabbit Polyclonal Antibody



CAB12792

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

## Applications:

WB

## Reactivity:

Arabidopsis

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

## Immunogen information

### Gene ID:

821625

### Uniprot

### Synonyms:

EIN3

### Immunogen:

A synthetic peptide of arabidopsis thaliana EIN3

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---