

CAB12813

Product Information

Product SKU:	CAB12813	Gene ID:	4892	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	197kDa	Conjugate:	Unconjugated
Calculated MW:	197kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable actin filament binding activity and muscle alpha-actinin binding activity. Predicted to be involved in cardiac muscle thin filament assembly. Predicted to be located in fascia adherens; muscle tendon junction; and myofibril. Predicted to be active in Z disc.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	N-RAP; NRAP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NRAP (NP_001248392.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.