

NRAP Rabbit Polyclonal Antibody



CAB12813

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

197kDa

Calculated MW:

192kDa/193kDa/197kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4892

Uniprot

Q86VF7

Synonyms:

NRAP; N-RAP


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NRAP (NP_001248392.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of SH-SY5Y cells, using NRAP antibody (CAB12813) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (CABM00021). Exposure time: 30s.