

CAB12833

Product Information

Product SKU:	CAB12833	Gene ID:	57488	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	102kDa	Conjugate:	Unconjugated
Calculated MW:	102kDa	Isotype:	IgG

Immunogen Information

Background:	Enables calcium ion binding activity; identical protein binding activity; and phospholipid binding activity. Predicted to be involved in endocytosis; endoplasmic reticulum-plasma membrane tethering; and lipid transport. Located in endoplasmic reticulum-plasma membrane contact site. Is extrinsic component of cytoplasmic side of plasma membrane; integral component of plasma membrane; and intrinsic component of endoplasmic reticulum membrane.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	E-Syt2; FAM62B; CHR2SYT
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 630-760 of human FAM62B (NP_065779.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.