

NKAIN4 Rabbit Polyclonal Antibody



CAB12841

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

23kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:1000 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

128414

Uniprot

Q8IVV8

Synonyms:

NKAIN4; C20orf58; FAM77A; bA261N11.2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 89-208 of human NKAIN4 (NP_690603.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
