

Product Information

Product SKU:	CAB12862	Gene ID:	84749	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	71kDa	Conjugate:	Unconjugated
Calculated MW:	59kDa	Isotype:	IgG

Immunogen Information

**Background:** USP30, a member of the ubiquitin-specific protease family (see USP1, MIM 603478), is a novel mitochondrial deubiquitinating (DUB) enzyme (Nakamura and Hirose, 2008 [PubMed 18287522]).

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
Synonyms:	USP30
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 57-390 of human USP30 (NP_116052.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.