

CAB12867

Product Information

Product SKU:	CAB12867	Gene ID:	94033	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	21kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ferric iron binding activity and ferrous iron binding activity. Involved in several processes, including cellular iron ion homeostasis; positive regulation of aconitate hydratase activity; and positive regulation of succinate dehydrogenase activity. Located in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MTF; FTMT
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 183-242 of human FTMT (NP_803431.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.