

CAB12875

Product Information

Product SKU:	CAB12875	Gene ID:	441478	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	12kDa	Conjugate:	Unconjugated
Calculated MW:	12kDa	Isotype:	IgG

Immunogen Information

Background: Involved in Notch signaling pathway. Acts upstream of or within negative regulation of Notch signaling pathway and positive regulation of canonical Wnt signaling pathway.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NRARP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human NRARP (NP_001004354.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.